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Conservser cette notice.

### en

## SHORT-STRETCH COMPRESSION BANDAGE

### Description/Destination

One-way short-stretch compression bandage, indicated in the treatment of vascular, lymphatic and traumatic conditions.

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the device dimensions.

**Composition**
Polyamide, elastan, viscosse, cotton.

### Properties/Mode of action

The bandages exert pressure on the limb, facilitating venous return by reducing venous dilation and edema volume.

### Indications

Treatment of edema (acute or chronic).

Treatment of lymphedema

Chronic venous insufficiency (as an alternative to compression stockings).

Following sclerotherapy or varicose vein surgery.

Treatment of venous ulcers (with ABI > 0.8 or under close medical supervision for an ABI entre 0.5 and 0.8).

Symptomatic treatment of deep or superficial acute venous thrombosis.

Prevention of post-thrombotic syndrome: alternative to compression stockings for at-risk patients;

Alternative to compression stockings for the treatment of edema of the legs or swelling during pregnancy.

Treatment of symptoms related to soft tissue injuries.

### Contraindications

Do not use in the event of known allergy to any of the components.

PAD (peripheral artery disease) with an ABI < 0.6.

Skin wounds, abscess, infected skin conditions of a limb concerned.

Advanced diabetic microangiopathy (if pressure > 30 mmHg)

Decompensated heart failure.

Phlegmasia coerulea dolens (painful blue phlebitis with arterial compression) of a limb concerned.

Extra-anatomic bypass.

Septic thrombosis of a limb concerned.

Severe peripheral neuropathy.

### Precautions

The choice of the width is at the discretion of the healthcare professional, depending on the limb to be treated.

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Strictly comply with your healthcare professional's prescription and recommendations for use.

Check the condition of the affected limb daily (with particular attention for patients with sensory deficit).

Do not apply the product directly to an open wound or apply a dressing in the event of damaged skin.

The bandage can be used with a foam bandage or a Varico in the event of dysmorphia or bony prominences.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

For hygiene and performance reasons, do not re-use the product for another patient.

### Undesirable side-effects

This device can cause skin reactions (redness, itching, burns, blisters, etc.) or wounds of various degrees of severity.

This device can cause pressure marks.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

### Instructions for use/Application

- Position the end of the bandage flat against the base of the toes.
- Wind round a first time going under the foot.
- The second wind should go towards the ankle. Cover the heel.
- Go under the foot and back under the ankle.
- Return to the top of the foot, pass underneath and go round the ankle.
- Go up the calf, spacing the turns regularly, following the prescription.

Check that there are no uncovered spaces between the turns.

### Care

Machine wash at 40°C, delicates cycle. Squeeze out excess water, if necessary. Can be autoclaved.

### Storage

Store at room temperature, preferably in the original packaging.

### Disposal

Dispose of in accordance with local regulations.

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